



Single-Patient Disposable Cuffs



Neonate Disposable Cuffs



Reusable Cuffs without Bladders



Reusable Cuffs with Bladders

- ◆ Comfortable high quality nylon material impregnated with polyurethane
- ◆ Latex-free, antibacterial coated and blood stain free
- ◆ Available in full range of AHA compatible size
- ◆ Meet AAMI standards for durability

NIBP Cuffs

Single-Patient Disposable Cuffs

Item	Part No.(Single Hose)	Part No.(Double Hoses)	Patient	Circumference	Length(mm)	Width(mm)
1	C1011	C1021	Infant	10-15 cm	194	57
2	C1012	C1022	Pediatric	14-21.5 cm	274	83
3	C1013	C1023	Small Adult	20.5-28cm	334	108
4	C1014	C1024	Adult	27-35cm	524	142
5	C1015	C1025	Adult-Long	27-35cm	605	142
6	C1016	C1026	Large Adult	34-43cm	608	170
7	C1017	C1027	Large Adult-Long	34-43cm	688	170
8	C1018	C1028	Thigh	42-54 cm	745	198

Neonate Disposable Cuffs

Item	Part No.(Single Hose)	Part No.(Double Hoses)	Patient	Circumference	Length(mm)	Width(mm)
1	C101A	C102A	Neonate	3.3-5.6cm	110	25
2	C101B	C102B	Neonate	4.2-7.1cm	140	32
3	C101C	C102C	Neonate	5.0-10.5cm	200	40
4	C101D	C102D	Neonate	6.9-11.7cm	250	50
5	C101E	C102E	Neonate	8.9-15cm	270	60

Multi-Patient Reusable Cuffs without Bladders

Item	Part No.(Single Hose)	Part No.(Double Hoses)	Patient	Circumference	Length(mm)	Width(mm)
1	C0011	C0021	Infant	10-15 cm	215	55
2	C0012	C0022	Pediatric	14-21.5 cm	270	75
3	C0013	C0023	Small Adult	20.5-28cm	410	115
4	C0014	C0024	Adult	27-35cm	520	140
5	C0015	C0025	Adult-Long	27-35cm	605	140
6	C0016	C0026	Large Adult	34-43cm	610	170
7	C0017	C0027	Large Adult-Long	34-43cm	688	170
8	C0018	C0028	Thigh	42-54 cm	770	215

Multi-Patient Reusable Cuffs with Bladders

Item	Part No.(Single Hose)	Part No.(Double Hoses)	Patient	Circumference	Length(mm)	Width(mm)
1	C0111	C0121	Infant	10-15 cm	215	65
2	C0112	C0122	Pediatric	14-21.5 cm	300	91
3	C0113	C0123	Small Adult	20.5-28cm	394	114
4	C0114	C0124	Adult	27-35cm	465	142
5	C0115	C0125	Adult-Long	27-35cm	600	142
6	C0116	C0126	Large Adult	34-43cm	553	174
7	C0117	C0127	Large Adult-Long	34-43cm	668	174
8	C0118	C0128	Thigh	42-54 cm	742	214

Solaris Medical Technology, Inc.



China Sales Office
 3/F, Zhongjian Building
 18 Yanshan Road, Nanshan District
 Shenzhen, Guangdong 518067 CHINA
 Tel: +86 755 26507042
 Fax: +86 755 26635815
 E-mail: sales.intl@szsolaris.com
 www.szsolaris.com

 3/F, Zhongjian Building
 18 Yanshan Road, Nanshan District
 Shenzhen, Guangdong 518067 CHINA
 Tel: +86 755 26635816
 Fax: +86 755 26635815
 E-mail: intl-sales@szsolaris.com
 www.szsolaris.com